

# PRODUCT CATALOGUE

.....Attaining Endless Possibilities

## About FATRA

FATRA was started in the late 2016 by a young technician, who had the skills to undertake the clinical maintenance of many different manufactured medical devices. Through innovation ideas and market driven economic savings, FATRA started to grow rapidly.

That progress continues today alongside the technological evolution and the development of clinical engineering, which could no longer be defined as just “biomedical engineering”, but as the advancement of clinical technology into an integrated management system within healthcare facilities.

01

### Our Vision

To increase the efficiency and improve the quality of services provided to citizens by health, social care facilities and public authorities.

02

### Our Mission

FATRA's mission is to develop and manage integrated outsourced clinical engineering and ICT solutions ensuring that all technologies used in hospital, health and social care Institutions..

03

### Our Aim

FATRA aims to offer innovative solutions to reduce healthcare spending in the technology sector and make this spending more focused.

## WHY CHOOSE US

Safety check report with manufacturer, model number, serial number and work performed.

We ensure detailed calibration/safety stickers on all tested equipment that we work on.

Our services ensure your equipment is delivering proper dosages and is safe for your

We ensure up to date and highest quality with our Computerized preventative maintenance

We take on patient care and physical medical equipment with utmost responsibility and

24 hour response time is typical and you can contact us anytime with troubleshooting questions

.....Attaining endless possibilities

## DIAMEDICA BABY CPAP 10



### CPAP 10

- Can be used for neonate up to 5 kg
- CPAP 10 twin flowmeters for up to 5l/min oxygen and 5l/min medical air
- Self-contained unit generates 95% oxygen and its own medical air
- Humidifier ensures appropriate warmth and moisture for the air and oxygen delivered to the patient.
- Clear oxygen/air mixing chart allows for easy setting of flow rates.
- Voltage protector is fitted to the CPAP machine, protecting the machine from unstable power surges

### CPAP 20

- Can be used for babies/pediatrics up to 20kg
- CPAP 10 twin flowmeters for up to 10l/min oxygen and 10l/min medical air
- Self-contained unit generates 95% oxygen and its own medical air
- Humidifier ensures appropriate warmth and moisture for the air and oxygen delivered to the patient.
- Clear oxygen/air mixing chart allows for easy setting of flow rates.
- Voltage protector is fitted to the CPAP machine, protecting the machine from unstable power surges

## DIAMEDICA PORTABLE BABY CPAP

### Features

- RAM cannula
- Flow rate of 8-20l/min dependent on cannula size
- up to 8cmH2O pressure dependent on nares seal.
- Supplementary oxygen port
- Passive humidification
- Battery life up to 8 hours at full use
- Recharge options: 12v vehicle or mains 110v-240v



03

.....Attaining endless possibilities

## SLE2000

Mechanical ventilator for infant.

### Features

- High performance constant flow
- Time Cycle Pressure Limited
- easy to set up
- easy to clean
- SLE valveless technology
- it allows for CMV, CPAP, PTV or SIMV
- High Frequency Oscillation.



## EKOM Dk50

Medical air compressor

### Features

- Motor performance: 0.43kW
- Air tank capacity: 4.7l
- Compressor efficiency: 50l/min/5bar
- Duty cycle: s1-100%
- Noise level: <53dB
- Noise reduction cabinet
- Weight: 31kg



04

.....Attaining endless possibilities

## DRIVE DeVilbiss

5Liters Oxygen concentrator



- Capacity: 5 liters
- Outlet Pressure: 8.5 psig
- Oxygen Percentage: 1-5 LPM, 87% - 96%
- Power Consumption: 310 Watt Average, 275 Watts @ 1.2 LPM and Below

## Drive DeVilbiss

10 Liter Oxygen Concentrator



- Capacity: 10 LPM
- Dimensions: 62.2x34.23x0.4 cm
- Sound level: Less than 40dB
- Outlet Pressure: 8.5 psig
- Oxygen Percentage: 1-5 LPM, 87% - 96%
- Power Consumption: 310 Watt Average, 275 Watts @ 1.2 LPM and Below

04

.....Attaining endless possibilities

## OMRON HEM-7121 STANDARD.

Automatic Blood Pressure Monitor

### Product Features

- Cuff Wrapping Guide. Guides user to the right cuff wrap by indicating "OK" to ensure an accurate and reliable blood pressure measurement.
- Stores 30 Measurements in Memory.
- Body Movement Detection. ...
- IntelliSense Technology. ...
- Irregular Heartbeat Detection. ...
- Hypertension Indicator.



## OMRON M2 CLASSIC.

Automatic Blood Pressure Monitor

### Product Features

- Cuff Wrapping Guide. Guides user to the right cuff wrap by indicating "OK" to ensure an accurate and reliable blood pressure measurement.
- Cuff Type: Medium Cuff 22-32 cm.
- Memory: Last Measurement
- Body Movement Detection.
- IntelliSense Technology.
- Irregular Heartbeat Detection.
- Hypertension Indicator.



05

.....Attaining endless possibilities



### OMRON M2 HEM-7121 BASIC. Automatic Blood Pressure Monitor

#### Product Features

- Cuff Wrapping Guide. Guides user to the right cuff wrap by indicating "OK" to ensure an accurate and reliable blood pressure measurement.
- Stores 30 Measurements in Memory.
- Body Movement Detection. ...
- IntelliSense Technology. ...
- Irregular Heartbeat Detection. ...
- Hypertension Indicator.



### OMRON 7 Series

Wireless Wrist blood pressure monitor

#### Product Features

- Cuff Wrapping Guide. Guides user to the right cuff wrap by indicating "OK" to ensure an accurate and reliable blood pressure measurement.
- Memory: 90 readings for one user or unlimited readings when used with the free OMRON Connect app
- The monitor fits adult wrists from 5.3" to 8.5" in circumference..
- Body Movement Detection.
- IntelliSense Technology.
- Irregular Heartbeat Detection..
- Hypertension Indicator.
- Unit operates on 2 AAA batteries (included)



.....Attaining endless possibilities

06

### BISTOS BT-720

Bistos Portable Patient Vital Signs Monitor.

#### FEATURES

- SpO2,
- Pulse,
- NIBP,
- The BT-720 has a 4.3" color touch screen.
- The PC software allows the vital signs information to be transferred and analyzed on a computer.



### BISTOS BT-710 4.3" . PULSE OXIMETER

#### FEATURES

- 4.3" color touch screen
- Handheld style
- SpO2, Pulse
- Accurate SpO2 performance during motion and low perfusion
- Specialized PC software for data review and analysis
- Over 8 hours continuous working on rechargeable lithium-ion battery
- 5V DC input with convenient Micro-USB charger
- Battery can be charged by external USB battery
- Suitable for adult, pediatric and neonate



.....Attaining endless possibilities

07

### BISOT BT-750

Multiparameter Patient's monitor

#### FEATURES

- Display
- ECG: 3/5 lead with printer
- Respiration
- SpO2
- CO2
- NIBP
- IBP
- Temperature
- Pulse rate
- Built-in battery: 6hours




### BISOT BT-780

Multiparameter Patient's monitor

#### FEATURES

- With a 15.6" LCD colour touch screen display
- Visual and acoustic alarm signals critical values
- Respiratory rate
- Non-invasive blood pressure (NIBP)
- Pulse rate
- Blood oxygen saturation (SpO2)
- Body temperature
- Physiological signals,
- ECG: 7 Lead
- With or without printer
- Rechargeable battery with 5 h operating time



.....Attaining endless possibilities

08

### ZONECARE PM-7000D

Multi-parameter Patient Monitor

#### FEATURES

- ECG
- NIBP
- SpO2
- Pulse rate
- Temperature
- EtCO2
- Built-in battery; 6hours
- Touch screen
- High resolution 12.1 inch screen design & intuitive shortcut keys
- Handle design convenient to carry
- Nurse call
- Night mode



### iMAC 1800

Professional 18-lead ECG machine

#### FEATURES

- 15-inch HD LCD Screen
- 18 Lead and 12 Lead printing format.
- Support pacemaker detection.
- Alphanumeric keyboard and tracking ball operation
- ECG data can be saved as PDF, JPG, PNG, XML etc
- Support data transmission to PC via LAN port, USB port, SD card and WiFi



.....Attaining endless possibilities

09

### ZONECARE iMAC 12

ECG MACHINE

#### FEATURES

- 12-lead synchronous acquisition and display
- 10.4-inch super sensitive and capacitive touch screen design
- A/D converter
- Support pacemaker detection
- Support data transmission to PC via LAN port, USB, SD card and WiFi
- Specialized storage bag
- Multifunction capacitive touch pen
- Built-in battery
- ECG data can be saved as HI7, XML, PDF, JPG and DICOM



### ZONECARE iMAC 300

Professional 3-Channel ECG machine

#### FEATURES

- 5-inch high resolution LCD screen
- 12-lead synchronous acquisition and display
- Support pacemaker detection.
- Alphanumeric keyboard and tracking ball operation
- ECG data can be saved as PDF, JPG, PNG, XML etc
- Support data transmission to PC via LAN port, USB port, SD card and WiFi
- Manual/Automatic printing.
- Support paperless records and storage electronic case.
- Stores more than 100 cases.



.....Attaining endless possibilities

### ZONECARE ViV 80

Ultrasound system

#### FEATURES

- Comprehensive clinical application
- B Steer, 3D/4D, Zpage, Zlive, PW, CW, AMM, CMM, TDI, TVI, PDI, DPDI.
- Panoramic imaging, Elastic imaging mode, EFOV, Biopsy function.
- Visualize user-defined presets function (Qsave)
- Ten user-defined programming keys.
- ECG port, USB 3.0, Built-in DVD
- DICOM 3.0
- Multiple data transmission
- SSD+HDD
- USB foot switch
- Built-in battery.
- Transducers: Convex, Linear, Transvaginal, Phased Array, Micro-convex, 4D Volume.



### ZONECARE ViV 60

Ultrasound system

#### FEATURES

- Comprehensive clinical application
- B Steer, 3D/4D, Zpage, Zlive, PW, CW, AMM, CMM, TDI, TVI, PDI, DPDI.
- Panoramic imaging, Elastic imaging mode, EFOV, Biopsy function.
- Visualize user-defined presets function (Qsave)
- Ten user-defined programming keys.
- ECG port, USB 3.0, Built-in DVD
- DICOM 3.0
- Multiple data transmission
- SSD+HDD
- USB foot switch
- Built-in battery.
- Transducers: Convex, Linear, Transvaginal, Phased Array, Micro-convex, 4D Volume.



.....Attaining endless possibilities

### ZONECARE U30

ECG MACHINE

#### FEATURES

- 12-lead synchronous acquisition and display
- 5-inch screen
- Digital sampling
- Support pacemaker detection
- Support data transmission to PC via LAN port, USB, SD card and WiFi
- Support Z-folded and roll thermal recording paper
- Built-in battery
- ECG data can be saved as HI7, XML, PDF, JPG and DICOM



### ZONECARE U90

Professional 12-Channel ECG machine

#### FEATURES

- 9.7-inch color touch screen and alphanumeric keyboard
- 12-lead synchronous acquisition and display
- Support pacemaker detection.
- Alphanumeric keyboard and tracking ball operation
- ECG data can be saved as PDF, JPG, PNG, XML etc
- Support data transmission to PC via LAN port, USB port, SD card and WiFi
- Manual/Automatic printing.
- Support paperless records and storage electronic case.
- Three - pronged lead wire, effectively reduce winding.



.....Attaining endless possibilities

### ZONECARE - N7

Color doppler Ultrasound system

#### FEATURES

- 19-inch high resolution medical LCD monitor
- Gel warmer
- Carotid IMT auto measurement.
- B Steer, 3D/4D, Zpage, Zlive, PW, CW, AMM, CMM, TDI, TVI, PDI, DPDI.
- Panoramic imaging, Elastic imaging mode, EFOV, Biopsy function.
- Visualize user-defined presets function (Qsave)
- Ten user-defined programming keys.
- ECG port, USB 3.0, Built-in DVD
- DICOM 3.0
- Multiple data transmission
- HDMI output
- USB foot switch
- Built-in battery.
- Transducers: Convex, Linear, Transvaginal, Phased Array, Micro-convex, 4D Volume.



### ZONECARE - Q3

Compact color doppler Ultrasound system

#### FEATURES

- 15-inch high resolution medical LCD monitor
- 10.4-inch touch screen.
- Gel warmer
- B Steer, 3D/4D, Zpage, Zlive, PW, CW, AMM, CMM, TDI, TVI, PDI, DPDI.
- Panoramic imaging, Elastic imaging mode, EFOV, Biopsy function.
- Visualize user-defined presets function (Qsave)
- Ten user-defined programming keys.
- ECG port, USB 3.0, Built-in DVD
- DICOM 3.0
- Multiple data transmission
- Optional USB foot switch
- USB foot switch
- Built-in battery.
- Transducers: Convex, Linear, Transvaginal, Phased Array, Micro-convex, 4D Volume.



.....Attaining endless possibilities

### ZONECARE - M5

Color doppler Ultrasound system

#### FEATURES

- 12-inch high resolution medical LCD monitor
- Two activated transducer connector.
- Carotid IMT auto measurement.
- B Steer, 3D/4D, Zpage, Zlive, PW, CW, HPRF, PDI, DPDI.
- Panoramic imaging, Elastic imaging mode, EFOV, Biopsy function.
- Visualize user-defined presets function (Qsave)
- Ten user-defined programming keys.
- ECG port, USB 3.0, Built-in DVD
- DICOM 3.0
- Multiple data transmission
- USB foot switch
- Built-in battery.
- Transducers: Convex, Linear, Transvaginal, Phased Array, Micro-convex.



### ZONECARE - i50

Portable B&W Ultrasound system

#### FEATURES

- 12-inch high resolution medical LCD monitor
- Two transducer connector.
- Auto freezing function
- Auto-generation BMP format report.
- G-guide
- One key data storage
- Full screen imaging FSI
- Multiple data transmission
- Optional USB foot switch
- USB foot switch
- Built-in battery.
- Transducers: Convex, Linear, Transvaginal, Phased Array, Micro-convex.



14

.....Attaining endless possibilities

### ZONECARE - iE75

Ambulatory blood pressure monitoring

#### FEATURES

- Support MS Windows 2003, Windows XP or higher version of Windows
- Support morning blood pressure index, variability index, arteriosclerosis index analysis.
- Analysis of various blood pressure parameters and statistical results.
- Automatically generate reports in PDF format
- Supports setting of various time interval measurement programs such as day, night, morning, special, etc
- Data transmission: Bluetooth or USB
- Adapt oscillometric method for blood pressure measurement.



### ZONECARE - IMAC 10 PC ECG

#### FEATURES

- Able to detect pacemaker
- Lead connection diagram
- Sampling rate is 32000HZ to make the waveform restoration more accurate.
- Bluetooth connection
- Resting ECG App (for Android)
  - Real Time Sampling
  - Automatic Diagnosis
  - Touch screen
- Workstation (for PC)
  - Printing
  - Waveform analysis
  - Average template.
  - Automatic Diagnosis
  - Measure
  - Reanalysis.



15

.....Attaining endless possibilities

### BABY INCUBATOR

Neonatal Intensive Care Unit Incubator

#### FEATURES

- Support Air temperature up to 37 degree celsius
- Support humidification
- Casters wheel for mobility.
- Two drop down doors for easier assess
- Support both manual and Servo control
- Tubing access ports: 13.
- Alarm system.



### GE GIRAFFE INCUBATOR

INFANT INCUBATOR

#### Features

- In-Bed scale.
- Ability to control heat, oxygen, light and noise.
- Patented rotating mattress to minimize baby stimulation
- Mattress size: 26 in x 19 in (66 cm x 48 cm)
- Tubing access ports: 13.
- Baby (servo) temperature control: 35-37.5° C in 0.1° increments.
- Air temperature control: 20-39° C in 0.1° increments.
- Servo humidity control range: 30-95% relative humidity in 5% increments.



16

.....Attaining endless possibilities

### ISOLLET C2000

Neonatal Intensive Care Unit Incubator

#### FEATURES

- Support Air temperature up to 37 degree celsius
- Support humidification
- Casters wheel for mobility.
- Two drop down doors for easier assess
- Support both manual and Servo control
- Tubing access ports: 13.
- Alarm system.
- Storage drawer.
- Drip stand.



### ISOLLET T1500

Dragger infant Transport Incubator

#### Features

- Quick access to the infant. The head access / slide out mattress feature allow quick and easy access to baby for intubation or CPR
- On casters, height-adjustable.
- With 4 operating windows.
- Weight. 49.2 kg, 72 kg (108.5 lb) ; Length. 95.9 cm, 102 cm (37.8 in).
- Double wall; Skin temperature probe; O2 inlet; Examination lamp; 2 access doors



17

.....Attaining endless possibilities





### CENTRIFUGE BUCKET CENTRIFUGE

#### Features:

- 1Capacity: 20ml\*6 ;
- Power source. 220V 50HZ 110V 60HZ ;
- Speed Range ( rpm). 0 ~ 4000 ( rpm ) ;
- Time Range(min). 0-60min.



### CENTRIFUGE BUCKET CENTRIFUGE

#### Features:

- 1Capacity: 20ml\*10 ;
- Power source. 220V 50HZ 110V 60HZ ;
- Speed Range ( rpm). 0 ~ 4000 ( rpm ) ;
- Time Range(min). 0-60min.

18

.....Attaining endless possibilities



### 80-2B HAEMATOCRIT CENTRIFUGE

#### Features:

- Rotor type: Disk rotor hematocrite rotor
- Speed max.: 13000 rpm
- Frequency (Centrifuges): 50/60 Hz
- Capacity of vessels: 0.2 ml
- Cooling (Centrifuges): Air cooling



### HC-12A HAEMATOCRIT CENTRIFUGE

#### Features:

- Rotor type: Disk rotor hematocrite rotor
- Speed max.: 13000 rpm
- Frequency (Centrifuges): 50/60 Hz
- Capacity of vessels: 0.2 ml
- Cooling (Centrifuges): Air cooling

19

.....Attaining endless possibilities



### ERMA EL-120 Electrolyte Analyzer

#### Features:

- Parameters: Na+, K+, Cl- Ca++
- Sample Type. Analyses whole blood, serum, plasma and urine
- Relative Humidity up to 85% non- condensing
- Power Supply: 100 to 240 VAC  $\pm$ 10%, 50Hz
- Dimensions: 395 x 260 x 375 mm
- Weight 7.5 kg Apprx

### M104 Electrolyte Analyzer.

#### Features.

- Compact, cost-efficient and easy - to- use.
- Excellent precision and reliability.
- Long life, high performance, maintenance -free electrodes.
- Automatic sampling, washing and calibration.
- Touch screen ,large LCD display.
- Internal thermal recorder ,
- supporting barcode reader.



20

.....Attaining endless possibilities



### Abbot i-STAT 1 Handheld Blood Gas Analyzer

#### Features:

- It has a rechargeable battery,
- can scan 2D barcodes,
- Offers on-screen guidance during the testing process.
- Identify problems and potential solutions in real time.

### ST-200CC ABGEM Blood Gas Analyzer

#### Features:

- The ST-200 CC ABG Analyzer measures various combinations of Sodium, Potassium, Ionized Calcium, Lithium, pH, Chloride, pCO<sub>2</sub>, pO<sub>2</sub>, Hct, and calculated parameters like Hb, HCO<sub>3</sub>, BE, BE-B, BE-ECF, TCO<sub>2</sub>, AG, AG(K), O<sub>2</sub> Sat, O<sub>2</sub> Ct and SBC in whole blood, serum, plasma, and urine.



21

.....Attaining endless possibilities



### Olympus Cx23

Clinical Light Microscope

- Features:**
- Objectives: 4x, 10x, 40x, 100x Oil
  - Condenser: 1.25 N.A.
  - Stage: Rackless mechanical stage
  - Focusing: Coaxial coarse/fine focus system with limit stopper and torque adjustment
  - Light: LED



### EDAN i15

Blood Gas Analyzer

- Features:**
- Calculated Values:  $cH^+$ ,  $HCO_3^-$ -act,  $HCO_3^-$ -std,  $BE(ecf)$ ,  $BE(B)$ ,  $BB(B)$ ,  $ctCO_2$ ,  $sO_2(est)$ ,  $Ca^{++}(7.4)$ ,  $AnGap$ ,  $tHb(est)$ ,
  - $pO_2(A-a)$ ,  $pO_2(a/A)$ ,  $RI$ ,  $pO_2/FIO_2$ ,  $cH^+(T)$ ,  $pH(T)$ ,  $pCO_2(T)$ ,  $pO_2(T)$ ,  $pO_2(A-a)(T)$ ,  $pO_2(a/A)(T)$ ,  $RI(T)$ ,  $pO_2(T)/FIO_2$ ,  $mOsm$



### IB-9025A

Laboratory Incubator

- Features:**
- Power input: 220V - 240V
  - Temperature: ...Celsius



### DY - 300

Electrophoresis Machine

- Features:**
- Genotype:

23

.....Attaining endless possibilities

22

.....Attaining endless possibilities



### POL - EKO

Laboratory Incubator

- Features:**
- Interior volume: up to 50 l 51 - 100 l 101 - 200 l 201 - 300 l 301 - 400 l 701 - 800 l 901 - 1000 l
  - Air convection forced
  - Load max.: 11 - 50 kg 51 - 100 kg 101 - 150 kg
  - Depth inside (incubators,...): 200 mm 250 mm 360 mm 450 mm 510 mm 560 mm 600 mm
  - Cabinet type: heating oven
  - Height inside: 230 mm 320 mm 390 mm 540 mm 720 mm 800 mm 1040 mm 1200 mm 1610 mm
  - Maximum temperature (heat): + 100 °C
  - Width inside: 320 mm 400 mm 460 mm 470 mm 600 mm 800 mm 1040 mm
  - Depth outside: 470 mm 520 mm 620 mm 710 mm 770 mm 780 mm 820 mm 870 mm 880 mm
  - Width outside: 510 mm 590 mm 650 mm 810 mm 1010 mm 1260 mm



### Reichert UNISTAT®

Bilirubinometer

- Features:**
- Measures TSB concentration in 5 seconds
  - Measures undiluted serum
  - Corrects for oxyhemoglobin
  - Full measuring range
  - Automatic calibration setting



### BM-100 Series

Transcutaneous Jaundice Detector

- Features:**
- 3" colour touch screen.
  - long battery life (2000 measurements from full charge)
  - storage of patient information and results.
  - long life span (> 150,000 scans)
  - self calibration.
  - bar code/ QR reader to scan baby ID/ nurse ID (BM-100C only)
  - USB transfer of data (BM-100C only).
  - Applicable to Neonatology, Pediatric, NICU etc



### BT-400

Neonatal Jaundice treatment.

- Features:**
- Light Source - Blue LEDs
  - Wavelength: Peak between 450 ~475 nm.
  - Intensity: 20 ~ 100W/CM<sup>2</sup>/nm at 40 cm.
  - Variation in intensity over 6 hours: ± 10% (within illumination area)
  - Effective surface area: 40 x 20 cm.

24

.....Attaining endless possibilities

25

.....Attaining endless possibilities



### NEOSoft NEUVIZ ACE Computer Tomography

#### Features:

- Unique 60kV low dose technology. offers optimizes dose level and image quality.
- Up-to-date spectral imaging.
- Application: for whole-body tomography
- Number of slices: 32-slice
- Real-time MPR
- Prism imaging.
- 2D, 3D, Virtual Endoscopy, MIP/AIP, Vessel analysis, Brain Perfusion, Body perfusion, Dental analysis, Lung nodules Analysis, Virtual Colonoscopy, Lung Density Evaluation, Fat Analysis, Tumor Assessment, Nerve system DSA, Lung Nodules ROI, and Bone Measurement.

### NEUVERSION 640 Floor Mount Digital X-RAY

#### Features:

- X-ray generator
- High sensitive anti-collision system.
- High precision real-time tracking system
- Patented image stitching technology
- Wireless Flat Panel Detector
- Full featured Imaging Software and Excellent Digital Image processing
- Provides convenient user interface and easy operation.
- Radiographic stand & automatic collimator control function.
- Low dose algorithm technology



### HB-7021 AUTOMATIC HEMATOLOGY ANALYZER

#### Features:

- Easy to use hematology analyzer.
- Open system with 60 test per hour performance.
- Defining 19 blood parameters with leukocytes differentiation on 3 subpopulations and building 3 histogram
- Able to operate with small-volume blood samples.
- Good choice for laboratories, who needs affordable analyzer with high performance.



### NEO-BIL PLUS BILIRUBINOMETER

#### Features:

- Measurement of bilirubin tests
- = Serum minimum volume: 10uL
- Reading cuvette: heparinized haematocrit capillary.
- Measure unit: mg/dL or  $\mu\text{mol/L}$
- Completely automated.
- Reading time: approximately 3 seconds.
- Result: on LCD and printer.
- Data storage: last 1000 result
- Interface: USB PORT

26

.....Attaining endless possibilities

27

.....Attaining endless possibilities



### INFANT RADIANT WARMER

#### Features:

- Both manual and Servo mode of heating
- LED illumination light.
- Caster wheel.
- Side trail.
- Heater up to 80 degree Celsius
- APGAR Score
- Temperature sensor

### DIGITAL RCS-20 ELECTRONIC INFANT WEIGHT SCALE

#### Features:

- Function: Body weight.
- Display type: LED
- Weight range: 1kg - 20kg
- Power: Battery



### WEIGHT AND HEIGHT SCALE

#### Features:

- Structure: Electronic Scale
- Display: Analog
- Measuring range: 70 - 190cm
- Maximum weight: 120kg
- Minimum scale: 0.5kg
- Minimum division: 0.5cm

### CLINICAL DENTAL UNIT

#### Features:

- LED oral light.
- Automatic water supply
- Stainless steel hand piece
- Stainless steel foot switch
- Safety switch ON/OFF
- Silent and stable motor
- Nine chair-positions memory.
- Double arm rest.
- Removable rotating spittoon.
- Display screen.



28

.....Attaining endless possibilities

29

.....Attaining endless possibilities



### PATIENT'S BED

#### Features:

- Caster
- Brake
- Side rail
- Three crank



### PATIENT'S BED ELECTRIC/MANUAL BED.

#### Features:

- Caster wheel
- Brake
- Siderail
- Bed frame
- Hand release
- Remote control.
- Four crank



### NANNY BABY APNEA MONITOR

#### Features:

- Easy to use
- Durable.
- Sensor pad up to 14cm deep





### AMBULANCE STRETCHER

#### Features:

- Caster wheel
- Siderail
- Bed frame
- Hand release
- Mattress
- Patient belt.




# INSTALLATIONS





# Installations done by our team



.....Attaining endless possibilities



# Installations done by our team



.....Attaining endless possibilities



# Installations done by our team



.....Attaining endless possibilities



# Installations done by our team



.....Attaining endless possibilities



# Repairs done by our team



.....Attaining endless possibilities




# Sourcing devices accessories



.....Attaining endless possibiities



## REPAIRS AND MAINTENANCE






## REPAIRS AND MAINTENANCE





## REPAIRS AND MAINTENANCE





# EQUIPMENT'S EXHIBITIONS






# TRAININGS





## SOME OF OUR CLIENTS



**Federal Medical  
Centre Asaba**  
Delta State, Nigeria



**THE  
OBIJACKSON**  
men & Children's Hospital




**LIFELINE**  
Children's Hospital

*.....Attaining endless possibilities*



**CREST CONSULTANT HOSPITAL  
AND FERTILITY CENTRE**



E-mail: amaitaribariweni@yahoo.

*.....Attaining endless possibilities*



**EDI INTERNATIONAL  
HOSPITAL LTD.**



**St. Philomena's  
Hospital**



**Our Mother Theresa Specialist  
Hospital & Fertility Center,  
Asaba**



**BENIN COOPERATIVE SOCIETY LTD**

*.....Attaining endless possibilities*



## Our Partners

---



---

**BIOBASE**

**zoncare**



.....Attaining endless possibilities



## CONTACTS:



102, Mission Road, Benin City, Edo State.



[fatramedical@gmail.com](mailto:fatramedical@gmail.com)



[fatramedicalequipment](https://www.instagram.com/fatramedicalequipment)



+2347054773403 +2349022057670



+2347054773403

.....Attaining endless possibilities